

AUDIOSCHOOL LITE

SIMPLIFIED AND INTUITIVE AUDIOMETRY



AUDIOMETRY

STAND-
ALONE



TOUCH
SCREEN



EASY
TO USE



RECHARGEABLE
BATTERY

This portable **audiometer** has been designed to allow you to perform **fast and accurate** audiometry tests in **ambulatory mode**. Its **touch screen** offers you the ability to carry out your measurements very easily.

Easy to use, it is intended for **school screening** and is suitable for both **doctors and school nurses**.

MADE IN FRANCE

TESTS PERFORMED ON THE DEVICE

Thanks to its carrying case, **Audioschool Lite** can easily be taken with you wherever you need it. Ultra light, it fits easily in your hand. Via its **touch screen**, you can **easily perform the tests** and then transfer the results to your audiogram notepad..



TECHNICAL CHARACTERISTICS

Pure tone audiometry measurements:

Headset: standard

Frequency range: 125 to 8 000 Hz

Sound intensity: -10 to 100 dB HL

Steps: 5 dB

Dual channel: yes

Stimulation indicator light

Type of stimulation: continuous or pulsed

Stimulation button option: keep holding or click

Device

Power supply: lithium battery (5 hours of autonomy), level displayed on the screen

Size: 9 x 11 x 3,6 cm

Weight: 239 g

Transport: bag (27 x 20 x 12 cm), less than 1 kg with the headset

Touch screen usable with finger or with the provided stylus



Audioschool Lite audiometer



Audiometry headset



USB cable and power adapter



Audiogram notepad



Carry-on bag

